

Human Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3) ELISA Kit | TRPV3 elisa kit

Human Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3) ELISA Kit

CATALOG #
AAA186941

REACTIVITY
Human


PRODUCT SYNONYMS

Transient Receptor Potential Cation Channel Subfamily V; Member 3 (TRPV3); Human Transient Receptor Potential Cation Channel Subfamily V; Member 3 (TRPV3) ELISA Kit; VRL3; Vanilloid receptor-like 3; TRPV3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3). No significant cross-reactivity or interference between Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3) and analogues was observed.
Sequence Length	765
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.781-50ng/mL
Sensitivity	< 0.29ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TRPV3 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Transient Receptor Potential Cation Channel Subfamily V, Member 3 (TRPV3) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	384551638
NCBI GeneID	162514
NCBI Accession #	NP_001245134.1
NCBI GenBank Nucleotide #	NM_001258205.1
UniProt Accession #	Q8NET8; Q8NDW7; Q8NET9; Q8NFH2
Molecular Weight	87,661 Da
NCBI Official Full Name	transient receptor potential cation channel subfamily V member 3 isoform 1
NCBI Official Synonym Full Names	transient receptor potential cation channel subfamily V member 3
NCBI Official Symbol	TRPV3
NCBI Official Synonym Symbols	OLMS; VRL3; FNEPPK2
NCBI Protein Information	transient receptor potential cation channel subfamily V member 3
UniProt Protein Name	Transient receptor potential cation channel subfamily V member 3
UniProt Gene Name	TRPV3
UniProt Synonym Gene Names	TrpV3; VRL-3

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.