

Mouse Periostin (POSTN) ELISA Kit | POSTN elisa kit

Mouse Periostin (POSTN) ELISA Kit

CATALOG #
AAA187297

REACTIVITY
Mouse


PRODUCT SYNONYMS

Periostin (POSTN); Mouse Periostin (POSTN) ELISA Kit; OSF-2; PDLPOSTN; PN; Osteoblast Specific Factor 2; POSTN elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of PER1. No significant cross-reactivity or interference between PER1 and analogues was observed.
Sequence Length	721
Samples	Human Tissue Homogenates, Cell Lysates And Other Biological Fluids
Assay Type	Quantitative Sandwich
Detection Range	15.6-1,000pg/mL
Sensitivity	5.5pg/mL.
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level PER1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level PER1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for POSTN elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of PER1 in human tissue homogenates, cell lysates and other biological fluids. Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to PER1. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to PER1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PER1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of PER1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	209862907
NCBI GeneID	10631
NCBI Accession #	NP_006466.2
NCBI GenBank Nucleotide #	NM_006475.2
UniProt Accession #	Q15063; Q15064; Q29XZ0; Q3KPJ5; Q5VSY5; Q8IZF9; B1ALD8; C0IMJ1; C0IMJ2; C0IMJ4; D2KRH7; F5H628
Molecular Weight	80,357 Da
NCBI Official Full Name	periostin isoform 1
NCBI Official Synonym Full Names	periostin
NCBI Official Symbol	POSTN
NCBI Official Synonym Symbols	PN; OSF2; OSF-2; PDLPOSTN
NCBI Protein Information	periostin
UniProt Protein Name	Periostin
UniProt Gene Name	POSTN
UniProt Synonym Gene Names	OSF2; PN; OSF-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.