

Human Guanine Deaminase (GDA) ELISA Kit | GDA elisa kit

Human Guanine Deaminase (GDA) ELISA Kit

CATALOG #
AAA187439

REACTIVITY
Human


PRODUCT SYNONYMS

Guanine Deaminase (GDA); Human Guanine Deaminase (GDA) ELISA Kit; CYPIN; Guanase; Guanine aminase; Guanine aminohydrolase; p51-nedasin; GDA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Guanine Deaminase (GDA). No significant cross-reactivity or interference between Guanine Deaminase (GDA) and analogues was observed.
Sequence Length	454
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue Homogenates and Other Biological Fluids
Detection Range	0.312-20ng/mL
Sensitivity	< 0.129ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Guanine Deaminase (GDA) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Guanine Deaminase (GDA) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GDA elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Guanine Deaminase (GDA). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Guanine Deaminase (GDA). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Guanine Deaminase (GDA), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Guanine Deaminase (GDA) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	148747414
NCBI GeneID	83585
NCBI Accession #	NP_113964.2
NCBI GenBank Nucleotide #	NM_031776.2
UniProt Accession #	Q9WTT6
Molecular Weight	51,016 Da
NCBI Official Full Name	guanine deaminase
NCBI Official Synonym Full Names	guanine deaminase
NCBI Official Symbol	Gda
NCBI Official Synonym Symbols	Cypin
NCBI Protein Information	guanine deaminase
UniProt Protein Name	Guanine deaminase
UniProt Gene Name	Gda
UniProt Synonym Gene Names	Guanase; Guanine aminase; GAH

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.