

Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit | CILP2 elisa kit

Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit

CATALOG #
AAA187527

REACTIVITY
Human


PRODUCT SYNONYMS

Cartilage Intermediate Layer Protein 2 (CILP2); Human Cartilage Intermediate Layer Protein 2 (CILP2) ELISA Kit; CILP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Cartilage Intermediate Layer Protein 2 (CILP2). No significant cross-reactivity or interference between Cartilage Intermediate Layer Protein 2 (CILP2) and analogues was observed.
Sequence Length	1,156
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.312-20ng/mL
Sensitivity	< 0.116ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Cartilage Intermediate Layer Protein 2 (CILP2) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Cartilage Intermediate Layer Protein 2 (CILP2) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CILP2 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cartilage Intermediate Layer Protein 2 (CILP2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Cartilage Intermediate Layer Protein 2 (CILP2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cartilage Intermediate Layer Protein 2 (CILP2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Cartilage Intermediate Layer Protein 2 (CILP2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	145699139
NCBI GeneID	148113
NCBI Accession #	NP_694953.2
NCBI GenBank Nucleotide #	NM_153221.2
UniProt Accession #	Q8IUL8; Q6NV88; Q8N4A6; Q8WV21
Molecular Weight	126,291 Da
NCBI Official Full Name	cartilage intermediate layer protein 2
NCBI Official Synonym Full Names	cartilage intermediate layer protein 2
NCBI Official Symbol	CILP2
NCBI Official Synonym Symbols	CLIP-2
NCBI Protein Information	cartilage intermediate layer protein 2
UniProt Protein Name	Cartilage intermediate layer protein 2
UniProt Gene Name	CILP2
UniProt Synonym Gene Names	CILP-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cartilage-intermediate-layer-protein-2-cilp2/187527>