

Human Coiled Coil Domain Containing Protein 80 (CCDC80) ELISA Kit | CCDC80 elisa kit

Human Coiled Coil Domain Containing Protein 80 (CCDC80) ELISA Kit

CATALOG #
AAA187610

REACTIVITY
Human


PRODUCT SYNONYMS

Coiled Coil Domain Containing Protein 80 (CCDC80); Human Coiled Coil Domain Containing Protein 80 (CCDC80) ELISA Kit; DRO1; SSG1; URB; Steroid Sensitive Gene 1; Down-regulated by oncogenes protein 1; Up-regulated in BRS-3 deficient mouse homolog; CCDC80 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Coiled Coil Domain Containing Protein 80 (CCDC80). No significant cross-reactivity or interference between Coiled Coil Domain Containing Protein 80 (CCDC80) and analogues was observed.
Sequence Length	85
Assay Type	Double-antibody Sandwich
Samples	Serum, Plasma, Tissue Homogenates and Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.054ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Coiled Coil Domain Containing Protein 80 (CCDC80) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Coiled Coil Domain Containing Protein 80 (CCDC80) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CCDC80 elisa kit	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Coiled Coil Domain Containing Protein 80 (CCDC80). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Coiled Coil Domain Containing Protein 80 (CCDC80). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Coiled Coil Domain Containing Protein 80 (CCDC80), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Coiled Coil Domain Containing Protein 80 (CCDC80) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	41152076
NCBI GeneID	151887
NCBI Accession #	NP_955806.1
NCBI GenBank Nucleotide #	NM_199512.2
UniProt Accession #	Q76M96; Q5PR20; Q6GPG9; Q8IVT6; Q8NBV1; Q8NHY8; D3DN67
Molecular Weight	10,174 Da
NCBI Official Full Name	coiled-coil domain-containing protein 80
NCBI Official Synonym Full Names	coiled-coil domain containing 80
NCBI Official Symbol	CCDC80
NCBI Official Synonym Symbols	URB; DRO1; SSG1; okuribin
NCBI Protein Information	coiled-coil domain-containing protein 80
UniProt Protein Name	Coiled-coil domain-containing protein 80
UniProt Gene Name	CCDC80
UniProt Synonym Gene Names	DRO1; URB

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/coiled-coil-domain-containing-protein-80-ccdc80/187610>