

Human Chaperonin Containing TCP1, Subunit 2 (CCT2) ELISA Kit | CCT2 elisa kit

Human Chaperonin Containing TCP1, Subunit 2 (CCT2) ELISA Kit

CATALOG #
AAA187702

REACTIVITY
Human


PRODUCT SYNONYMS

Chaperonin Containing TCP1; Subunit 2 (CCT2); Human Chaperonin Containing TCP1; Subunit 2 (CCT2) ELISA Kit; CCT-beta; CCTB; TCP-1-beta; T-Complex Protein 1 Subunit Beta; CCT2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Chaperonin Containing TCP1, Subunit 2 (CCT2). No significant cross-reactivity or interference between Chaperonin Containing TCP1, Subunit 2 (CCT2) and analogues was observed.
Sequence Length	493
Assay Type	Double-antibody Sandwich
Samples	Tissue Homogenates, Cell Lysates and Other Biological Fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.056ng/mL
Application	Enzyme-linked immunosorbent assay for Antigen Detection.
Intra-assay Precision (Precision within an assay)	3 samples with low, middle and high level Chaperonin Containing TCP1, Subunit 2 (CCT2) were tested 20 times on one plate, respectively.
Inter-assay Precision (Precision between assays)	3 samples with low, middle and high level Chaperonin Containing TCP1, Subunit 2 (CCT2) were tested on 3 different plates, 8 replicates in each plate.
CV(%) =	SD/meanX100
Intra-Assay	CV<10%
Inter-Assay	CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CCT2 elisa kit

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Chaperonin Containing TCP1, Subunit 2 (CCT2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Chaperonin Containing TCP1, Subunit 2 (CCT2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Chaperonin Containing TCP1, Subunit 2 (CCT2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Chaperonin Containing TCP1, Subunit 2 (CCT2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	119600607
NCBI GeneID	10693
UniProt Accession #	Q92526; Q8TC34; B4DX20; B4DYB0
Molecular Weight	53,586 Da
NCBI Official Full Name	chaperonin containing TCP1, subunit 6B (zeta 2)
NCBI Official Synonym Full Names	chaperonin containing TCP1 subunit 6B
NCBI Official Symbol	CCT6B
NCBI Official Synonym Symbols	Cctz2; CCTZ-2; TSA303; CCT-zeta-2; TCP-1-zeta-2
NCBI Protein Information	T-complex protein 1 subunit zeta-2
UniProt Protein Name	T-complex protein 1 subunit zeta-2
UniProt Gene Name	CCT6B
UniProt Synonym Gene Names	TCP-1-zeta-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/chaperonin-containing-tcp1-subunit-2-cct2/187702>