

Human Nitric Oxide Synthase 2, Inducible (NOS2) ELISA Kit | NOS2 elisa kit

Human Nitric Oxide Synthase 2, Inducible (NOS2) High Sensitive ELISA Kit

CATALOG # AAA188523	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Nitric Oxide Synthase 2; Inducible (NOS2); Human Nitric Oxide Synthase 2; Inducible (NOS2) High Sensitive ELISA Kit; Inducible; NOS2A; INOS; HEP-NOS; I-NOS; Hepatocytes Oxide Synthase; Peptidyl-cysteine S-nitrosylase NOS2; NOS2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2). No significant cross-reactivity or interference between High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2) and analogues was observed.
Sequence Length	1153
Assay Type	Double-antibody Sandwich
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Detection Range	15.6-1,000pg/mL
Sensitivity	The minimum detectable dose of this kit is typically less than 5.5pg/mL
Applicable Applications for NOS2 elisa kit	ELISA
Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2) were tested 20 times on one plate, respectively. Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%
Test Principle	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of High Sensitive Nitric Oxide Synthase 2, Inducible (NOS2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	24041029
NCBI GeneID	4843
NCBI Accession #	NP_000616.3
NCBI GenBank Nucleotide #	NP_000616.3
UniProt Accession #	P35228
NCBI Official Full Name	nitric oxide synthase, inducible
NCBI Official Synonym Full Names	nitric oxide synthase 2
NCBI Official Symbol	NOS2
NCBI Official Synonym Symbols	NOS; INOS; NOS2A; HEP-NOS
NCBI Protein Information	nitric oxide synthase, inducible
UniProt Protein Name	Nitric oxide synthase, inducible
UniProt Gene Name	NOS2
UniProt Synonym Gene Names	NOS2A; HEP-NOS; Inducible NOS; iNOS
UniProt Entry Name	NOS2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/nitric-oxide-synthase-2-inducible-nos2/188523>