

Human Inhibitory Subunit Of NF Kappa B Alpha (IκBa) ELISA Kit | IκBa elisa kit

Human Inhibitory Subunit Of NF Kappa B Alpha (IκBa) High Sensitive ELISA Kit

CATALOG # AAA188556	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Inhibitory Subunit Of NF Kappa B Alpha (IκBa); Human Inhibitory Subunit Of NF Kappa B Alpha (IκBa) High Sensitive ELISA Kit; Inhibitory Subunit Of NF Kappa B Alpha; NFKBIA; IKB-A; MAD-3; NFKBI; Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-Cells Inhibitor; Alpha; Major histocompatibility complex enhancer-binding protein MAD3; IκBa elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa). No significant cross-reactivity or interference between High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa) and analogues was observed.
Sequence Length	317
Assay Type	Double-antibody Sandwich
Samples	Tissue homogenates, cell lysates and other biological fluids
Detection Range	31.2-2,000pg/mL
Sensitivity	The minimum detectable dose of this kit is typically less than 11.4pg/mL
Applicable Applications for IκBa elisa kit	ELISA
Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa) were tested 20 times on one plate, respectively. Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%
Test Principle	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of High Sensitive Inhibitory Subunit Of NF Kappa B Alpha (IκBa) in the samples is then determined by comparing the O.D. of the samples to the standard curve.
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	10092619
NCBI GeneID	4792
NCBI Accession #	NP_065390.1
NCBI GenBank Nucleotide #	NP_065390.1
UniProt Accession #	P25963
NCBI Official Full Name	NF-kappa-B inhibitor alpha
NCBI Official Synonym Full Names	NFKB inhibitor alpha
NCBI Official Symbol	NFKBIA
NCBI Official Synonym Symbols	IKBA; MAD-3; NFKBI; EDAID2
NCBI Protein Information	NF-kappa-B inhibitor alpha
UniProt Protein Name	NF-kappa-B inhibitor alpha
UniProt Gene Name	NFKBIA
UniProt Synonym Gene Names	IKBA; MAD3; NFKBI; Ikb-alpha
UniProt Entry Name	IKBA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/inhibitory-subunit-of-nf-kappa-b-alpha-ikba/188556>