

Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) ELISA Kit | NTXI elisa kit

Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) Mini Samples ELISA Kit


CATALOG # AAA188710	REACTIVITY Mouse	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------

PRODUCT SYNONYMS

Cross Linked N-Telopeptide Of Type I Collagen (NTXI); Mouse Cross Linked N-Telopeptide Of Type I Collagen (NTXI) Mini Samples ELISA Kit; Cross Linked N-Telopeptide Of Type I Collagen; NTX-I; NTX1; NTXI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI). No significant cross-reactivity or interference between Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) and analogues was observed.
Assay Type	Double-antibody Sandwich
Samples	Serum, plasma, urine and other biological fluids
Detection Range	2.47-200ng/mL
Sensitivity	The minimum detectable dose of this kit is typically less than 0.89ng/mL
Applicable Applications for NTXI elisa kit	ELISA
Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) were tested 20 times on one plate, respectively. Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%
Test Principle	This assay employs the competitive inhibition enzyme immunoassay technique. A monoclonal antibody specific to Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) has been pre-coated onto a microplate. A competitive inhibition reaction is launched between biotin labeled Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) and unlabeled Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) (Standards or samples) with the pre-coated antibody specific to Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI). After incubation the unbound conjugate is washed off. Next, avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. The amount of bound HRP conjugate is reverse proportional to the concentration of Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) in the sample. After addition of the substrate solution, the intensity of color developed is reverse proportional to the concentration of Mini Samples Cross Linked N-Telopeptide Of Type I Collagen (NTXI) in the sample.
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.



Standard Curve (Sample)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-elisa-kits/cross-linked-n-telopeptide-of-type-i-collagen-ntxi/188710>