

Human N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) CLIA Kit | NT-ProBNP clia kit

Human N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) CLIA Kit

CATALOG # AAA188919	REACTIVITY Human	APPLICATIONS CLIA (Chemiluminescence Immunoassay)
-------------------------------	----------------------------	---


PRODUCT SYNONYMS

N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP); Human N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) CLIA Kit; N-Terminal Pro-Brain Natriuretic Peptide; NT-Pro-BNP; ; N-BNP; NT-ProBNP clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP). No significant cross-reactivity or interference between N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) and analogues was observed.
Sequence Length	134
Assay Type	Double-antibody Sandwich
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Detection Range	3.43-2,500pg/mL
Sensitivity	The minimum detectable dose of this kit is typically less than 1.33pg/mL
Applicable Applications for NT-ProBNP clia kit	CLIA (Chemiluminescence Immunoassay)
Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) were tested 20 times on one plate, respectively. Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%
Test Principle	The microplate provided in this kit has been pre-coated with an antibody specific to N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the N-Terminal Pro-Brain Natriuretic Peptide (NT-ProBNP) level in the sample or standard.
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	4505433
NCBI GeneID	4879
NCBI Accession #	NP_002512.1
NCBI GenBank Nucleotide #	NP_002512.1
UniProt Accession #	P16860
NCBI Official Full Name	natriuretic peptides B preproprotein
NCBI Official Synonym Full Names	natriuretic peptide B
NCBI Official Symbol	NPPB
NCBI Official Synonym Symbols	BNP
NCBI Protein Information	natriuretic peptides B
UniProt Protein Name	Natriuretic peptides B
UniProt Gene Name	NPPB
UniProt Synonym Gene Names	BNP(1-32); BNP-32
UniProt Entry Name	ANFB_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/n-terminal-pro-brain-natriuretic-peptide-nt-probnp/188919>