

## Human Tissue Factor (TF) CLIA Kit | TF clia kit

Human Tissue Factor (TF) CLIA Kit

|                               |                            |                                                             |
|-------------------------------|----------------------------|-------------------------------------------------------------|
| CATALOG #<br><b>AAA188923</b> | REACTIVITY<br><b>Human</b> | APPLICATIONS<br><b>CLIA (Chemiluminescence Immunoassay)</b> |
|-------------------------------|----------------------------|-------------------------------------------------------------|


### PRODUCT SYNONYMS

Tissue Factor (TF); Human Tissue Factor (TF) CLIA Kit; Tissue Factor; CD142; FIII; F3; TFA; Thromboplastin; Coagulation Factor III; TF clia kit

### PRODUCT SPECIFICATIONS

|                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                              | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Specificity                             | This assay has high sensitivity and excellent specificity for detection of Tissue Factor (TF). No significant cross-reactivity or interference between Tissue Factor (TF) and analogues was observed.                                                                                                                                                                                                                                                                                                                                                                                     |
| Sequence Length                         | 238                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Assay Type                              | Double-antibody Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Samples                                 | plasma                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Detection Range                         | 2.74-2,000pg/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Sensitivity                             | The minimum detectable dose of this kit is typically less than 0.97pg/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Applicable Applications for TF clia kit | CLIA (Chemiluminescence Immunoassay)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Precision                               | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Tissue Factor (TF) were tested 20 times on one plate, respectively.<br>Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Tissue Factor (TF) were tested on 3 different plates, 8 replicates in each plate.<br>CV(%) = SD/meanX100<br>Intra-Assay: CV<10%<br>Inter-Assay: CV<12%                                                                                                                                                                    |
| Test Principle                          | The microplate provided in this kit has been pre-coated with an antibody specific to Tissue Factor (TF). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Tissue Factor (TF). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Tissue Factor (TF) level in the sample or standard. |
| Preparation and Storage                 | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.                                                                                                |

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

|                                  |                                       |
|----------------------------------|---------------------------------------|
| NCBI GI #                        | 296010912                             |
| NCBI GeneID                      | 2152                                  |
| NCBI Accession #                 | NP_001171567.1                        |
| NCBI GenBank Nucleotide #        | NP_001171567.1                        |
| UniProt Accession #              | P13726                                |
| NCBI Official Full Name          | tissue factor isoform 2               |
| NCBI Official Synonym Full Names | coagulation factor III, tissue factor |
| NCBI Official Symbol             | F3                                    |
| NCBI Official Synonym Symbols    | TF; TFA; CD142                        |
| NCBI Protein Information         | tissue factor                         |
| UniProt Protein Name             | Tissue factor                         |
| UniProt Gene Name                | F3                                    |
| UniProt Synonym Gene Names       | TF                                    |
| UniProt Entry Name               | TF_HUMAN                              |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/tissue-factor-tf/188923>