

# Human Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) CLIA Kit | VEGFR1 clia kit

Human Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) CLIA Kit

|                               |                            |                                                             |
|-------------------------------|----------------------------|-------------------------------------------------------------|
| CATALOG #<br><b>AAA189044</b> | REACTIVITY<br><b>Human</b> | APPLICATIONS<br><b>CLIA (Chemiluminescence Immunoassay)</b> |
|-------------------------------|----------------------------|-------------------------------------------------------------|


## PRODUCT SYNONYMS

Vascular Endothelial Growth Factor Receptor 1 (VEGFR1); Human Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) CLIA Kit; Vascular Endothelial Growth Factor Receptor 1; Flt1; Fms-Related Tyrosine Kinase 1; Vascular Permeability Factor Receptor; Fms-like tyrosine kinase 1; Tyrosine-protein kinase receptor FLT; VEGFR1 clia kit

## PRODUCT SPECIFICATIONS

|                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|----------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reactivity</b>                                  | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Specificity</b>                                 | This assay has high sensitivity and excellent specificity for detection of Vascular Endothelial Growth Factor Receptor 1 (VEGFR1).<br>No significant cross-reactivity or interference between Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) and analogues was observed.                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Sequence Length</b>                             | 687                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Assay Type</b>                                  | Double-antibody Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Samples</b>                                     | Serum, plasma, tissue homogenates, saliva, cell lysates, follicular fluid and other biological fluids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Detection Range</b>                             | 27.4-20,000pg/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Sensitivity</b>                                 | The minimum detectable dose of this kit is typically less than 10.5pg/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Applicable Applications for VEGFR1 clia kit</b> | CLIA (Chemiluminescence Immunoassay)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Precision</b>                                   | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) were tested 20 times on one plate, respectively.<br>Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) were tested on 3 different plates, 8 replicates in each plate.<br>CV(%) = SD/meanX100<br>Intra-Assay: CV<10%<br>Inter-Assay: CV<12%                                                                                                                                                                                                        |
| <b>Test Principle</b>                              | The microplate provided in this kit has been pre-coated with an antibody specific to Vascular Endothelial Growth Factor Receptor 1 (VEGFR1). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Vascular Endothelial Growth Factor Receptor 1 (VEGFR1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Vascular Endothelial Growth Factor Receptor 1 (VEGFR1) level in the sample or standard. |
| <b>Preparation and Storage</b>                     | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.                                                                                                                                                                                                            |

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

|                                  |                                                         |
|----------------------------------|---------------------------------------------------------|
| NCBI GI #                        | 229892220                                               |
| NCBI GeneID                      | 2321                                                    |
| NCBI Accession #                 | NP_001153392.1                                          |
| NCBI GenBank Nucleotide #        | NP_001153392.1                                          |
| UniProt Accession #              | P17948                                                  |
| NCBI Official Full Name          | vascular endothelial growth factor receptor 1 isoform 2 |
| NCBI Official Synonym Full Names | fms related receptor tyrosine kinase 1                  |
| NCBI Official Symbol             | FLT1                                                    |
| NCBI Official Synonym Symbols    | FLT; FLT-1; VEGFR1; VEGFR-1                             |
| NCBI Protein Information         | vascular endothelial growth factor receptor 1           |
| UniProt Protein Name             | Vascular endothelial growth factor receptor 1           |
| UniProt Gene Name                | FLT1                                                    |
| UniProt Synonym Gene Names       | FLT; FRT; VEGFR1; VEGFR-1; FLT-1; FLT                   |
| UniProt Entry Name               | VGFR1_HUMAN                                             |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.