

## Mouse Perilipin 4 (PLIN4) CLIA Kit | PLIN4 clia kit

Mouse Perilipin 4 (PLIN4) CLIA Kit

|                               |                            |                                                             |
|-------------------------------|----------------------------|-------------------------------------------------------------|
| CATALOG #<br><b>AAA189103</b> | REACTIVITY<br><b>Mouse</b> | APPLICATIONS<br><b>CLIA (Chemiluminescence Immunoassay)</b> |
|-------------------------------|----------------------------|-------------------------------------------------------------|


### PRODUCT SYNONYMS

Perilipin 4 (PLIN4); Mouse Perilipin 4 (PLIN4) CLIA Kit; Perilipin 4; Adipocyte protein S3-12; PLIN4 clia kit

### PRODUCT SPECIFICATIONS

|                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|--------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                                 | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Specificity                                | This assay has high sensitivity and excellent specificity for detection of Perilipin 4 (PLIN4). No significant cross-reactivity or interference between Perilipin 4 (PLIN4) and analogues was observed.                                                                                                                                                                                                                                                                                                                                                                                      |
| Assay Type                                 | Double-antibody Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Samples                                    | Tissue homogenates, cell lysates and other biological fluids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Detection Range                            | 1.37-1,000pg/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Sensitivity                                | The minimum detectable dose of this kit is typically less than 0.55pg/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Applicable Applications for PLIN4 clia kit | CLIA (Chemiluminescence Immunoassay)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Precision                                  | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Perilipin 4 (PLIN4) were tested 20 times on one plate, respectively.<br>Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Perilipin 4 (PLIN4) were tested on 3 different plates, 8 replicates in each plate.<br>$CV(\%) = SD/mean \times 100$<br>Intra-Assay: $CV < 10\%$<br>Inter-Assay: $CV < 12\%$                                                                                                                                                 |
| Test Principle                             | The microplate provided in this kit has been pre-coated with an antibody specific to Perilipin 4 (PLIN4). Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to Perilipin 4 (PLIN4). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the Perilipin 4 (PLIN4) level in the sample or standard. |
| Preparation and Storage                    | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.                                                                                                   |

### IMAGES



Standard Curve (Sample)

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/mouse-clia-kits/perilipin-4-plin4/189103>