

Human Chorionic Gonadotropin Beta Polypeptide (CGb) ELISA Kit | CGb elisa kit

Anti-Chorionic Gonadotropin Beta Polypeptide Antibody (Anti-CGb) ELISA Kit

CATALOG # AAA189238	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Chorionic Gonadotropin Beta Polypeptide (CGb); Anti-Chorionic Gonadotropin Beta Polypeptide Antibody (Anti-CGb) ELISA Kit; CG-B; CGB3; hCGB; HCG; CGb5; CGb7; CGb8; Choriogonadotropin subunit beta; CGb elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Anti-CGb. No significant cross-reactivity or interference between Anti-CGb and analogues was observed.
Samples	Serum, plasma and other biological fluids
Detection Range	3.12-200ng/mL
Sensitivity	< 1.22ng/mL
Applicable Applications for CGb elisa kit	ELISA
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Anti-CGb were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Anti-CGb were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CGB elisa kit	Intended Uses: The kit is an enzyme immunoassay for in vitro quantitative measurement of Anti-CGb in human serum, plasma and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with CGb. Standards or samples are then added to the appropriate microtiter plate wells with a Horseradish Peroxidase (HRP)-conjugated secondary antibody. After TMB substrate solution is added, those wells that contain Anti-CGb will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of Anti-CGb in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	1027923599
NCBI GeneID	1082
UniProt Accession #	P0DN86; P01233; Q13991; Q14000; Q3KPI3; Q3SY41; Q8WTT5; Q8WXL1; Q8WXL2; Q8WXL3; A1A5E0; B9ZVP5
Molecular Weight	19,053 Da
NCBI Official Full Name	Choriogonadotropin subunit beta 3
NCBI Official Synonym Full Names	chorionic gonadotropin beta subunit 3
NCBI Official Symbol	CGB3
NCBI Official Synonym Symbols	CGB; CGB5; CGB7; CGB8; hCGB
NCBI Protein Information	choriogonadotropin subunit beta 3
UniProt Protein Name	Choriogonadotropin subunit beta 3
UniProt Gene Name	CGB3
UniProt Synonym Gene Names	CGB; CG-beta

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/chorionic-gonadotropin-beta-polypeptide-cgb/189238>