

## Human Vav 1 Oncogene (VAV1) ELISA Kit | VAV1 elisa kit

Vav 1 Oncogene (VAV1) ELISA Kit

|                               |                            |                              |
|-------------------------------|----------------------------|------------------------------|
| CATALOG #<br><b>AAA189375</b> | REACTIVITY<br><b>Human</b> | APPLICATIONS<br><b>ELISA</b> |
|-------------------------------|----------------------------|------------------------------|


### PRODUCT SYNONYMS

Vav 1 Oncogene (VAV1); Vav 1 Oncogene (VAV1) ELISA Kit; Proto Oncogene Vav; Vav 1 Guanine Nucleotide Exchange Factor; VAV1 elisa kit

### PRODUCT SPECIFICATIONS

|                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                                 | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity                                | This assay has high sensitivity and excellent specificity for detection of VAV1. No significant cross-reactivity or interference between VAV1 and analogues was observed.                                                                                                                                                                                                                                                                                                                  |
| Samples                                    | Human tissue homogenates, cell lysates and other biological fluids                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Assay Type                                 | Quantitative Sandwich                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Detection Range                            | 0.156-10ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Sensitivity                                | < 0.057ng/mL                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Applicable Applications for VAV1 elisa kit | ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Intra-assay Precision                      | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level VAV1 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%                                                                                                                                                                                                                                                                                                                     |
| Inter-assay Precision                      | Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level VAV1 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%                                                                                                                                                                                                                                                                                   |
| Preparation and Storage                    | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.<br>To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. |

### IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

|                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Related Product Information for VAV1 elisa kit</b> | <p>Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of VAV1 in human tissue homogenates, cell lysates and other biological fluids.</p> <p>Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to VAV1. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to VAV1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain VAV1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of VAV1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.</p> |
|-------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/vav-1-oncogene-vav1/189375>