

Human Ubiquitin Like Modifier Activating Enzyme 2 (UBA2) ELISA Kit | UBA2 elisa kit

Human Ubiquitin Like Modifier Activating Enzyme 2 (UBA2)(Sandwich ELISA) ELISA Kit

CATALOG # AAA189586	REACTIVITY Human	APPLICATIONS ELISA
-------------------------------	----------------------------	------------------------------


PRODUCT SYNONYMS

Ubiquitin Like Modifier Activating Enzyme 2 (UBA2); Human Ubiquitin Like Modifier Activating Enzyme 2 (UBA2)(Sandwich ELISA) ELISA Kit; ELISA Kit for Ubiquitin Like Modifier Activating Enzyme 2 (UBA2); UBA2/Ubiquitin Like Modifier Activating Enzyme 2; SAE2; ARX; SUMO1 Activating Enzyme Subunit 2; Anthracycline-associated resistance ARX; Ubiquitin-like 1-activating enzyme E1B; SUMO-activating enzyme subunit 2; UBA2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of UBA2. No significant cross-reactivity or interference between UBA2 and analogues was observed.
Assay Type	Quantitative Sandwich
Samples	Human tissue homogenates, cell lysates and other biological fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.059ng/mL
Applicable Applications for UBA2 elisa kit	ELISA
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level UBA2 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level UBA2 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Store at 4 degree C for 1 month/ -20 degree C for long-term (within shelf life). The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for UBA2 elisa kit	Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to UBA2. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to UBA2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain UBA2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of UBA2 in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of UBA2 in human tissue homogenates, cell lysates and other biological fluids.
Product Categories/Family for UBA2 elisa kit	Endocrinology

NCBI AND UNIPROT INFORMATION

NCBI GI #	4885649
NCBI GeneID	10054
NCBI Accession #	NP_005490.1
NCBI GenBank Nucleotide #	NM_005499.2
UniProt Accession #	Q9UBT2; O95605; Q59H87; Q6IBP6; Q9NTJ1; Q9UED2
Molecular Weight	71,224 Da
NCBI Official Full Name	SUMO-activating enzyme subunit 2
NCBI Official Synonym Full Names	ubiquitin-like modifier activating enzyme 2
NCBI Official Symbol	UBA2
NCBI Official Synonym Symbols	ARX; SAE2; HRIHFB2115
NCBI Protein Information	SUMO-activating enzyme subunit 2; SUMO1 activating enzyme subunit 2; SUMO-1 activating enzyme subunit 2; ubiquitin-like 1-activating enzyme E1B; anthracycline-associated resistance ARX; UBA2, ubiquitin-activating enzyme E1 homolog
UniProt Protein Name	SUMO-activating enzyme subunit 2
UniProt Gene Name	UBA2
UniProt Synonym Gene Names	SAE2; UBLE1B
UniProt Entry Name	SAE2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/ubiquitin-like-modifier-activating-enzyme-2-uba2/189586>