

Human Complement Component 4d (C4d) CLIA Kit | C4d clia kit

Human CLIA Kit for Complement Component 4d (C4d) CLIA Kit

CATALOG #
AAA189719

REACTIVITY
Human


PRODUCT SYNONYMS

Complement Component 4d (C4d); Human CLIA Kit for Complement Component 4d (C4d) CLIA Kit; C4d/Complement Component 4d; C4d clia kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of C4d. No significant cross-reactivity or interference between C4d and analogues was observed.
Samples	Human serum, plasma and other biological fluids
Detection Range	1.37-1,000pg/mL
Sensitivity	< 0.61pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level C4d were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level C4d were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Store at 4 degree C for 1 month/-20 degree C for long-term (within shelf life) Shelf Life: 8 months The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for C4d clia kit	Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to C4d. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to C4d. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then the mixture of substrate A and B is added to generate glow light emission kinetics. Upon plate development, the intensity of the emitted light is proportional to the C4d level in the sample or standard. Intended Uses: The Chemiluminescent Immunoassay kit is designed for the in vitro sensitive quantitative measurement of C4d in human serum, plasma and other biological fluids.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	178557739
NCBI GeneID	721
NCBI Accession #	NP_001002029.3
NCBI GenBank Nucleotide #	NM_001002029.3
UniProt Accession #	P0C0L5; P01028; P78445; Q13160; Q13906; Q14033; Q14835; Q6U2E9; Q6U2G1; Q6U2I5; Q6U2L1; A2BHY4
Molecular Weight	192,751 Da
NCBI Official Full Name	complement C4-B preproprotein
NCBI Official Synonym Full Names	complement component 4B (Chido blood group)
NCBI Official Symbol	C4B
NCBI Official Synonym Symbols	CH; C4F; CO4; C4B1; C4B2; C4B3; C4B5; C4BD; C4B12; C4B_2; CPAMD3
NCBI Protein Information	complement C4-B; C3 and PZP-like alpha-2-macroglobulin domain-containing protein 3; Chido form of C4; basic complement C4; complement C4B1a
UniProt Protein Name	Complement C4-B
UniProt Gene Name	C4B
UniProt Synonym Gene Names	CO4; CPAMD3
UniProt Entry Name	CO4B_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-clia-kits/complement-component-4d-c4d/189719>