

Human NADPH Oxidase 4 ELISA Kit | NOX4 elisa kit

Human NADPH Oxidase 4 (NOX4) (Sandwich ELISA) Mini sample ELISA Kit

CATALOG #
AAA189953

REACTIVITY
Human


PRODUCT SYNONYMS

NADPH Oxidase 4; Human NADPH Oxidase 4 (NOX4) (Sandwich ELISA) Mini sample ELISA Kit; KOX; KOX-1; RENOX; Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4; Mini Samples ELISA Kit for NADPH Oxidase 4 (NOX4); NOX4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of NOX4. No significant cross-reactivity or interference between NOX4 and analogues was observed.
Assay Type	Quantitative Sandwich
Samples	Human tissue homogenates, cell lysates and other biological fluids
Detection Range	0.156-10ng/mL
Sensitivity	< 0.055ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level NOX4 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level NOX4 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	4 degree C for 1 month/ -20 degree C for long-term(within shelf life). Shelf life is one year

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NOX4 elisa kit	Principle of the Assay: The microplate provided in this kit has been pre-coated with an antibody specific to NOX4. Standards or samples are then added to the appropriate microplate wells with a biotin-conjugated antibody specific to NOX4. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain NOX4, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of NOX4 in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of NOX4 in human tissue homogenates, cell lysates and other biological fluids.
Product Categories/Family for NOX4 elisa kit	Signal transduction

NCBI AND UNIPROT INFORMATION

NCBI GI #	219842348
NCBI GeneID	50507
NCBI Accession #	NP_001137309.1
NCBI GenBank Nucleotide #	NM_001143837.1
UniProt Accession #	Q9NPH5; Q5K3R4; Q5K3R5; Q5K3R6; Q5K3R8; Q7Z7G3; Q86V92; A8K715
Molecular Weight	66,932 Da
NCBI Official Full Name	NADPH oxidase 4 isoform c
NCBI Official Synonym Full Names	NADPH oxidase 4
NCBI Official Symbol	NOX4
NCBI Official Synonym Symbols	KOX; KOX-1; RENOX
NCBI Protein Information	NADPH oxidase 4; kidney oxidase-1; renal NAD(P)H-oxidase; kidney superoxide-producing NADPH oxidase
UniProt Protein Name	NADPH oxidase 4
UniProt Gene Name	NOX4
UniProt Synonym Gene Names	RENOX; KOX-1
UniProt Entry Name	NOX4_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.