

Human Gasdermin D, GSDMD ELISA Kit | GSDMD elisa kit

Human Gasdermin D, GSDMD ELISA Kit

CATALOG #
AAA19096

REACTIVITY
Human


PRODUCT SYNONYMS

Gasdermin D; GSDMD; Human Gasdermin D; GSDMD ELISA Kit; GSDMD elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	484
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids
Detection Range	0.05-15ng/ml
Sensitivity	0.017ng/ml
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GSDMD elisa kit

Principle of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human GSDMD antibody. GSDMD present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human GSDMD Antibody is added and binds to GSDMD in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated GSDMD antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human GSDMD. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

Intended Uses: This sandwich kit is for the accurate quantitative detection of Human Gasdermin D (also known as GSDMD) in serum, plasma, cell culture supernates, Ascites, tissue homogenates or other biological fluids.

NCBI AND UNIPROT INFORMATION

NCBI GI #	261878494
NCBI GeneID	79792
NCBI GenBank Nucleotide #	NP_001159709.1
UniProt Accession #	P57764
NCBI Official Full Name	gasdermin-D
NCBI Official Synonym Full Names	gasdermin D
NCBI Official Symbol	GSDMD
NCBI Official Synonym Symbols	DF5L; DFNA5L; FKSG10; GSDMDC1
NCBI Protein Information	gasdermin-D
UniProt Protein Name	Gasdermin-D
UniProt Gene Name	GSDMD
UniProt Synonym Gene Names	DFNA5L; GSDMDC1
UniProt Entry Name	GSDMD_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/gasdermin-d-gsdmd/19096>