

Mouse Oxidized low-density lipoprotein ELISA Kit | OxLDL elisa kit

Mouse Oxidized low-density lipoprotein ELISA Kit

CATALOG #
AAA192089

REACTIVITY
Mouse


PRODUCT SYNONYMS

Oxidized low-density lipoprotein; Mouse Oxidized low-density lipoprotein ELISA Kit; OxLDL; OxLDL elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	Recombinant and natural Mouse Oxidized low-density lipoprotein
Samples	Serum, plasma, tissue homogenates, cell culture supernates or other biological fluids
Detection Range	78-5000 pg/mL
Intra-assay Precision	Intra-Assay CV: <=4.6%
Inter-assay Precision	Inter-Assay CV: <=6.9%
Preparation and Storage	For long term storage, please store the entire kit at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for OxLDL elisa kit

Intended Uses: This immunoassay kit allows for the in vitro quantitative determination of target antigen concentrations in serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.

Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to target antigen. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific for target antigen and then avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then a TMB substrate solution is added to each well. Only those wells that contain target antigen, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The concentration of target antigen in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/oxldl-mouse-elisa-kits/oxidized-low-density-lipoprotein/192089>