

Human NF-kappa-B essential modulator ELISA Kit | IKBKG elisa kit

Human NF-kappa-B essential modulator ELISA Kit

CATALOG #
AAA192648

REACTIVITY
Human


PRODUCT SYNONYMS

NF-kappa-B essential modulator; Human NF-kappa-B essential modulator ELISA Kit; NEMO; FIP-3; Ikb kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKBKG elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	487
Assay Type	Sandwich
Detection Range	0.78-50 ng/mL
Sensitivity	0.34ng/mL
Intra-Assay CV	<=4.8%
Inter-Assay CV	<=8.6%
Recovery	105%
Preparation and Storage	For long term storage, please store the entire kit at -20 degree C.

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	223671870
NCBI GeneID	8517
NCBI Accession #	NP_001093326.2
NCBI GenBank Nucleotide #	NM_001099856.3
UniProt Accession #	Q9Y6K9; Q7LBY6; Q7Z7F1
Molecular Weight	36,953 Da

NCBI Official Full Name	NF-kappa-B essential modulator isoform b
NCBI Official Synonym Full Names	inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase gamma
NCBI Official Symbol	IKBKG
NCBI Official Synonym Symbols	IP; IP1; IP2; FIP3; IPD2; NEMO; FIP-3; Fip3p; IMD33; AMCBX1; ZC2HC9; IKK-gamma
NCBI Protein Information	NF-kappa-B essential modulator
UniProt Protein Name	NF-kappa-B essential modulator
UniProt Gene Name	IKBKG
UniProt Synonym Gene Names	FIP3; NEMO; NEMO; IKKAP1; I-kappa-B kinase subunit gamma; IKK-gamma; IKKG; Ikb kinase subunit gamma
UniProt Entry Name	NEMO_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/ikbkg-human-elisa-kits/nf-kappa-b-essential-modulator/192648>