

Human Advanced glycation end-product ELISA Kit | AGE elisa kit

Human Advanced glycation end-product ELISA Kit

CATALOG #
AAA194553

REACTIVITY
Human


PRODUCT SYNONYMS

Advanced glycation end-product; Human Advanced glycation end-product ELISA Kit; AGE; AGE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	Recombinant and natural Human Advanced glycation end-product
Samples	Serum, plasma, tissue homogenates, cell culture supernates or other biological fluids
Detection Range	0.78-50 ng/mL
Sensitivity	< 0.41 ng/mL
Intra-assay Precision	<=5.3%
Inter-assay Precision	<=8.7%
Preparation and Storage	For long term storage, please store the entire kit at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AGE elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to target antigen. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific for target antigen and then avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. Then a TMB substrate solution is added to each well. Only those wells that contain target antigen, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of a sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The concentration of target antigen in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: This immunoassay kit allows for the in vitro quantitative determination of target antigen concentrations in serum, plasma, tissue homogenates, cell culture supernates or other biological fluids.
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/age-human-elisa-kits/advanced-glycation-end-product/194553>