

## Mouse Nuclear factor of activated T-cells 5 ELISA Kit | Nfat5 elisa kit

Mouse Nuclear factor of activated T-cells 5 ELISA Kit

|                               |                            |
|-------------------------------|----------------------------|
| CATALOG #<br><b>AAA196465</b> | REACTIVITY<br><b>Mouse</b> |
|-------------------------------|----------------------------|


### PRODUCT SYNONYMS

Nuclear factor of activated T-cells 5; Mouse Nuclear factor of activated T-cells 5 ELISA Kit; Nfat5 elisa kit

### PRODUCT SPECIFICATIONS

|                 |               |
|-----------------|---------------|
| Reactivity      | Mouse         |
| Sequence Length | 1415          |
| Assay Type      | Sandwich      |
| Detection Range | 78-5000 pg/ml |
| Sensitivity     | 40 pg/mL      |
| Intra-Assay CV  | <=6.6%        |
| Inter-Assay CV  | <=10.3%       |
| Recovery        | 105%          |

### IMAGES



Standard Curve (Sample)

### NCBI AND UNIPROT INFORMATION

|                                  |                                        |
|----------------------------------|----------------------------------------|
| NCBI GI #                        | <a href="#">41054033</a>               |
| NCBI GeneID                      | <a href="#">334390</a>                 |
| NCBI Accession #                 | <a href="#">NP_956189.1</a>            |
| NCBI GenBank Nucleotide #        | <a href="#">NM_199895.1</a>            |
| NCBI Official Full Name          | nuclear factor of activated T-cells 5  |
| NCBI Official Synonym Full Names | nuclear factor of activated T cells 5a |
| NCBI Official Symbol             | nfat5a                                 |
| NCBI Official Synonym Symbols    | nfat5; zgc:63502; wu:fi46f11           |
| NCBI Protein Information         | nuclear factor of activated T-cells 5  |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

