

Rabbit anti-General Sialic Acid (SA) Polyclonal Antibody | anti-SA antibody

Biotin-Linked Antibody to Sialic Acid (SA)

CATALOG # AAA20088	HOST Rabbit	REACTIVITY General	APPLICATIONS ELISA, Immunohistochemistry, Immunocytochemistry, Immunofluorescence
------------------------------	-----------------------	------------------------------	---


PRODUCT SYNONYMS

Sialic Acid (SA); Polyclonal Antibody; Biotin-Linked Antibody to Sialic Acid (SA); anti-SA antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	General
Clonality	Polyclonal
Purity/Purification	Antigen-specific Affinity Chromatography.
Form/Format	Supplied as solution form in 0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.
Concentration	0.58mg/mL (varies by lot)
Applicable Applications for anti-SA antibody	ELISA, CLIA, IHC (Immunohistochemistry), ICC (Immunocytochemistry), IF (Immunofluorescence)
Application Notes	Immunohistochemistry: 5-20ug/mL; Immunocytochemistry: 5-20ug/mL; Optimal working dilutions must be determined by end user.
Source	Antibody labeling
Label	Biotin
Traits	Liquid
Immunogen	Small Molecule, SA conjugated to OVA .
Preparation and Storage	Storage: Avoid repeated freeze/thaw cycles. Store at 4 degree C for frequent use. Aliquot and store at -20 degree C for 24 months. Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37 degree C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

IMAGES


WB (Western Blot)

Western Blot;Sample: Recombinant TNFa, Caprine.

ADDITIONAL INFORMATION

Product Categories/Family for anti-SA antibody	Biotin Antibody
--	-----------------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/polyclonal/sialic-acid-sa-antibody-immunohistochemistry-immunocytochemistry-immunofluorescence-elisa/20088>