

Human 25-Dihydroxyvitamin D3 (25(OH)D3/HVD3) ELISA Kit | 25(OH)D3/HVD3 elisa kit

Human 25-Dihydroxyvitamin D3 (25(OH)D3/HVD3) ELISA Kit

CATALOG #
AAA201984

REACTIVITY
Human

PRODUCT SYNONYMS

25-Dihydroxyvitamin D3 (25(OH)D3/HVD3); Human 25-Dihydroxyvitamin D3 (25(OH)D3/HVD3) ELISA Kit; 25(OH)D3/HVD3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	3ng/ml-48ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for 25(OH)D3/HVD3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of 25(OH)D3;HVD3 in the sample, This 25(OH)D3;HVD3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus 25(OH)D3;HVD3 concentration. The concentration of 25(OH)D3;HVD3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/25-dihydroxyvitamin-d3-25-oh-d3-hvd3/201984>