

Human C-Jun N-Terminal Kinases (JNK) ELISA Kit | JNK elisa kit

Human C-Jun N-Terminal Kinases (JNK) ELISA Kit

CATALOG #
AAA202062

REACTIVITY
Human


PRODUCT SYNONYMS

C-Jun N-Terminal Kinases (JNK); Human C-Jun N-Terminal Kinases (JNK) ELISA Kit; JNK elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	413
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1.5ng/ml-24ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for JNK elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of JNK in the sample, This JNK ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus JNK concentration. The concentration of JNK in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	325297084
NCBI GenBank Nucleotide #	NP_001191547.1
NCBI Official Full Name	c-Jun N-terminal kinase JNK

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/c-jun-n-terminal-kinases-jnk/202062>