

Human Melanoma Differentiation-Associated Gene 5 (MDA-5) ELISA Kit | MDA-5 elisa kit

Human Melanoma Differentiation-Associated Gene 5 (MDA-5) ELISA Kit

CATALOG #
AAA202585

REACTIVITY
Human

PRODUCT SYNONYMS

Melanoma Differentiation-Associated Gene 5 (MDA-5); Human Melanoma Differentiation-Associated Gene 5 (MDA-5) ELISA Kit; MDA-5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1002
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MDA-5 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration ofMDA-5in the sample, This MDA-5ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versusMDA-5concentration. The concentration ofMDA-5in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	304434770
NCBI GenBank Nucleotide #	NP_001182108.1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/melanoma-differentiation-associated-gene-5-mda-5/202585>