

Human Matrix Metalloproteinase 2 (MMP-2) ELISA Kit | MMP-2 elisa kit

Human Matrix Metalloproteinase 2 (MMP-2) ELISA Kit

CATALOG #
AAA202633

REACTIVITY
Human


PRODUCT SYNONYMS

Matrix Metalloproteinase 2 (MMP-2); Human Matrix Metalloproteinase 2 (MMP-2) ELISA Kit; MMP-2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	610
Sensitivity	0.1 ng/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MMP-2 elisa kit

Intended Uses: This MMP-2 ELISA kit is intended Laboratory for Research use only and is not for use in test or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MMP-2 in the sample, this MMP-2 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MMP-2 concentration. The concentration of MMP-2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	189217853
NCBI GeneID	4313
NCBI Accession #	NP_001121363.1

NCBI GenBank Nucleotide #	NP_001121363.1
UniProt Accession #	P08253
NCBI Official Full Name	72 kDa type IV collagenase isoform 2
NCBI Official Synonym Full Names	matrix metalloproteinase 2
NCBI Official Symbol	MMP2
NCBI Official Synonym Symbols	CLG4; MONA; CLG4A; MMP-2; TBE-1; MMP-II
NCBI Protein Information	72 kDa type IV collagenase
UniProt Protein Name	72 kDa type IV collagenase
UniProt Gene Name	MMP2
UniProt Synonym Gene Names	CLG4A; MMP-2
UniProt Entry Name	MMP2_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/matrix-metalloproteinase-2-mmp-2/202633>