

Human Tissue Inhibitors Of Metalloproteinase 3 (TIMP-3) ELISA Kit | TIMP-3 elisa kit

Human Tissue Inhibitors Of Metalloproteinase 3 (TIMP-3) ELISA Kit

CATALOG #
AAA202641

REACTIVITY
Human


PRODUCT SYNONYMS

Tissue Inhibitors Of Metalloproteinase 3 (TIMP-3); Human Tissue Inhibitors Of Metalloproteinase 3 (TIMP-3) ELISA Kit; TIMP-3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	31
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	10ng/ml-160ng/ml
Sensitivity	1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TIMP-3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TIMP-3 in the sample, This TIMP-3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TIMP-3 concentration. The concentration of TIMP-3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	16923457
NCBI Official Full Name	tissue inhibitor of metalloproteinase 3, partial

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tissue-inhibitors-of-metalloproteinase-3-timp-3/202641>