

Human A Disintegrin And Metalloprotease 5 (ADAMTS5) ELISA Kit | ADAMTS5 elisa kit

Human A Disintegrin And Metalloprotease 5 (ADAMTS5) ELISA Kit

CATALOG #
AAA202706

REACTIVITY
Human


PRODUCT SYNONYMS

A Disintegrin And Metalloprotease 5 (ADAMTS5); Human A Disintegrin And Metalloprotease 5 (ADAMTS5) ELISA Kit; ADAMTS5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	930
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	6.25ng/ml-100ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADAMTS5 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of ADAMTS5 in the sample, This ADAMTS5 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus ADAMTS5 concentration. The concentration of ADAMTS5 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	195539372
NCBI GeneID	11096
NCBI GenBank Nucleotide #	NP_008969.2
UniProt Accession #	Q9UNA0
NCBI Official Full Name	A disintegrin and metalloproteinase with thrombospondin motifs 5 preproprotein
NCBI Official Synonym Full Names	ADAM metalloproteinase with thrombospondin type 1 motif 5
NCBI Official Symbol	ADAMTS5
NCBI Official Synonym Symbols	ADMP-2; ADAM-TS5; ADAMTS-5; ADAMTS11; ADAM-TS 5; ADAMTS-11; ADAM-TS 11

NCBI Protein Information	A disintegrin and metalloproteinase with thrombospondin motifs 5
UniProt Protein Name	A disintegrin and metalloproteinase with thrombospondin motifs 5
UniProt Gene Name	ADAMTS5
UniProt Synonym Gene Names	ADAMTS11; ADMP2; ADAM-TS 5; ADAM-TS5; ADAMTS-5; ADAM-TS 11; ADAMTS-11
UniProt Entry Name	ATS5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/a-disintegrin-and-metalloprotease-5-adamts5/202706>