

Human Urokinase-Type Plasminogen Activator Receptor (uPAR) ELISA Kit | uPAR elisa kit

Human Urokinase-Type Plasminogen Activator Receptor (uPAR) ELISA Kit

CATALOG #
AAA202975

REACTIVITY
Human


PRODUCT SYNONYMS

Urokinase-Type Plasminogen Activator Receptor (uPAR); Human Urokinase-Type Plasminogen Activator Receptor (uPAR) ELISA Kit; uPAR elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	335
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for uPAR elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of uPAR in the sample, This uPAR ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus uPAR concentration. The concentration of uPAR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4505865
NCBI GeneID	5329
NCBI Accession #	NP_002650.1

NCBI GenBank Nucleotide #	NP_002650.1
UniProt Accession #	Q03405
NCBI Official Full Name	urokinase plasminogen activator surface receptor isoform 1
NCBI Official Synonym Full Names	plasminogen activator, urokinase receptor
NCBI Official Symbol	PLAUR
NCBI Official Synonym Symbols	CD87; UPAR; URKR; U-PAR
NCBI Protein Information	urokinase plasminogen activator surface receptor
UniProt Protein Name	Urokinase plasminogen activator surface receptor
UniProt Gene Name	PLAUR
UniProt Synonym Gene Names	MO3; UPAR; U-PAR; uPAR
UniProt Entry Name	UPAR_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/urokinase-type-plasminogen-activator-receptor-upar/202975>