

Human Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit | BPI elisa kit

Human Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit

CATALOG #
AAA203089

REACTIVITY
Human


PRODUCT SYNONYMS

Bactericidal/Permeability Increasing Protein (BPI); Human Bactericidal/Permeability Increasing Protein (BPI) ELISA Kit; BPI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	486
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BPI elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of BPI in the sample, This BPI ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus BPI concentration. The concentration of BPI in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	29244434
NCBI GeneID	329547
NCBI Accession #	NP_808518.1

NCBI GenBank Nucleotide #	NP_808518.1
UniProt Accession #	Q67E05
NCBI Official Full Name	bactericidal permeability-increasing protein isoform 1
NCBI Official Synonym Full Names	bactericidal permeability increasing protein
NCBI Official Symbol	Bpi
NCBI Official Synonym Symbols	Bpifd1; 9230105K17Rik
NCBI Protein Information	bactericidal permeability-increasing protein
UniProt Protein Name	Bactericidal permeability-increasing protein
UniProt Gene Name	Bpi
UniProt Synonym Gene Names	BPI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/bactericidal-permeability-increasing-protein-bpi/203089>