

Human Mitogen Activated Protein Kinase 3 (MAPK3) ELISA Kit | MAPK3 elisa kit

Human Mitogen Activated Protein Kinase 3 (MAPK3) ELISA Kit

CATALOG #
AAA203163

REACTIVITY
Human

PRODUCT SYNONYMS

Mitogen Activated Protein Kinase 3 (MAPK3); Human Mitogen Activated Protein Kinase 3 (MAPK3) ELISA Kit; MAPK3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	380
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	10 pg/ml-160 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MAPK3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of MAPK3 in the sample, This MAPK3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus MAPK3 concentration. The concentration of MAPK3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	8393331
NCBI GeneID	50689

NCBI GenBank Nucleotide #	NP_059043.1
UniProt Accession #	P21708
NCBI Official Full Name	mitogen-activated protein kinase 3
NCBI Official Synonym Full Names	mitogen activated protein kinase 3
NCBI Official Symbol	Mapk3
NCBI Official Synonym Symbols	p44; ERK1; ERT2; MNK1; Erk-1; Esrk1; MAPK1; Prkm3; p44erk1; p44mapk
NCBI Protein Information	mitogen-activated protein kinase 3
UniProt Protein Name	Mitogen-activated protein kinase 3
UniProt Gene Name	Mapk3
UniProt Synonym Gene Names	Erk1; Prkm3; MAP kinase 3; MAPK 3; ERK-1; p44-MAPK

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/mitogen-activated-protein-kinase-3-mapk3/203163>