

Human Triggering Receptor Expresses On Myeloid Cells-1 (TREM-1) ELISA Kit | TREM-1 elisa kit

Human Triggering Receptor Expresses On Myeloid Cells-1 (TREM-1) ELISA Kit

CATALOG #
AAA203178

REACTIVITY
Human


PRODUCT SYNONYMS

Triggering Receptor Expresses On Myeloid Cells-1 (TREM-1); Human Triggering Receptor Expresses On Myeloid Cells-1 (TREM-1) ELISA Kit; TREM-1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	225
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	15 pg/ml-240 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TREM-1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TREM-1 in the sample, This TREM-1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TREM-1 concentration. The concentration of TREM-1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

336088680

NCBI GeneID	54210
NCBI Accession #	NP_001229518.1
NCBI GenBank Nucleotide #	NP_001229518.1
UniProt Accession #	Q9NP99
NCBI Official Full Name	triggering receptor expressed on myeloid cells 1 isoform 2
NCBI Official Synonym Full Names	triggering receptor expressed on myeloid cells 1
NCBI Official Symbol	TREM1
NCBI Official Synonym Symbols	CD354; TREM-1
NCBI Protein Information	triggering receptor expressed on myeloid cells 1
UniProt Protein Name	Triggering receptor expressed on myeloid cells 1
UniProt Gene Name	TREM1
UniProt Synonym Gene Names	TREM-1
UniProt Entry Name	TREM1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/triggering-receptor-expresses-on-myeloid-cells-1-trem-1/203178>