

## Human Cytotoxin-Associated Protein (CagA) ELISA Kit | CagA elisa kit

Human Cytotoxin-Associated Protein (CagA) ELISA Kit

CATALOG #  
**AAA203272**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Cytotoxin-Associated Protein (CagA); Human Cytotoxin-Associated Protein (CagA) ELISA Kit; CagA elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Human                                                             |
| Sequence Length | 411                                                               |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 125 pg/ml-2000 pg/ml                                              |
| Sensitivity     | 10 pg/ml                                                          |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for CagA elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CAGA in the sample, This CAGAE LISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CAGA concentration. The concentration of CAGA in the samples is then determined by comparing the O.D. of the samples to the standard curve.

### NCBI AND UNIPROT INFORMATION

|                         |               |
|-------------------------|---------------|
| NCBI GI #               | 22335896      |
| NCBI Official Full Name | CagA, partial |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cytotoxin-associated-protein-caga/203272>