

Human Signal-Induced Proliferation-Associated Protein-1 (SPA-1/SIPA-1) ELISA Kit | SPA-1/SIPA-1 elisa kit

Human Signal-Induced Proliferation-Associated Protein-1 (SPA-1/SIPA-1) ELISA Kit

CATALOG #
AAA203339

REACTIVITY
Human


PRODUCT SYNONYMS

Signal-Induced Proliferation-Associated Protein-1 (SPA-1/SIPA-1); Human Signal-Induced Proliferation-Associated Protein-1 (SPA-1/SIPA-1) ELISA Kit; SPA-1/SIPA-1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1042
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	100 pg/ml-1600 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SPA-1/SIPA-1 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SPA-1;SIPA-1 in the sample, This SPA-1;SIPA-1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SPA-1;SIPA-1 concentration. The concentration of SPA-1;SIPA-1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

206725526

NCBI GeneID	6494
NCBI GenBank Nucleotide #	NP_694985.29
UniProt Accession #	Q96FS4
NCBI Official Full Name	signal-induced proliferation-associated protein 1
NCBI Official Synonym Full Names	signal-induced proliferation-associated 1
NCBI Official Symbol	SIPA1
NCBI Official Synonym Symbols	SPA1
NCBI Protein Information	signal-induced proliferation-associated protein 1
UniProt Protein Name	Signal-induced proliferation-associated protein 1
UniProt Gene Name	SIPA1
UniProt Synonym Gene Names	SPA1; Sipa-1
UniProt Entry Name	SIPA1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/signal-induced-proliferation-associated-protein-1-spa-1-sipa-1/203339>