

Human TARC/CCL17 ELISA Kit | TARC/CCL17 elisa kit

Human TARC/CCL17 ELISA Kit

CATALOG #
AAA203345

REACTIVITY
Human


PRODUCT SYNONYMS

TARC/CCL17; Human TARC/CCL17 ELISA Kit; TARC/CCL17 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	94
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	50 pg/ml-800 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TARC/CCL17 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TARC;CCL17 in the sample, This TARC;CCL17 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TARC;CCL17 concentration. The concentration of TARC;CCL17 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	3024711
NCBI GeneID	6361

UniProt Accession #	Q92583
NCBI Official Full Name	C-C motif chemokine 17
NCBI Official Synonym Full Names	C-C motif chemokine ligand 17
NCBI Official Symbol	CCL17
NCBI Official Synonym Symbols	TARC; ABCD-2; SCYA17; A-152E5.3
NCBI Protein Information	C-C motif chemokine 17
UniProt Protein Name	C-C motif chemokine 17
UniProt Gene Name	CCL17
UniProt Synonym Gene Names	SCYA17; TARC
UniProt Entry Name	CCL17_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tarc-ccl17/203345>