

Human Endothelium-Derived Hyperpolarizing Factor (EDHF) ELISA Kit | EDHF elisa kit

Human Endothelium-Derived Hyperpolarizing Factor (EDHF) ELISA Kit

CATALOG #
AAA203363

REACTIVITY
Human


PRODUCT SYNONYMS

Endothelium-Derived Hyperpolarizing Factor (EDHF); Human Endothelium-Derived Hyperpolarizing Factor (EDHF) ELISA Kit; EDHF elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sensitivity	1 pg/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for EDHF elisa kit

Intended Uses: This EDHF ELISA kit is intended Laboratory for Research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of EDHF in the sample, this EDHF ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus EDHF concentration. The concentration of EDHF in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/endothelium-derived-hyperpolarizing-factor-edhf/203363>