

Human Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 (Tie-2) ELISA Kit | Tie-2 elisa kit

Human Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 (Tie-2) ELISA Kit

CATALOG #
AAA203377

REACTIVITY
Human


PRODUCT SYNONYMS

Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 (Tie-2); Human Tyrosine Kinase With Immunoglobulin-Like And EGF-Like Domains 2 (Tie-2) ELISA Kit; Tie-2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	150 pg/ml-2400 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Tie-2 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Tie-2 in the sample, This Tie-2 Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Tie-2 concentration. The concentration of Tie-2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	121946267
UniProt Accession #	Q8IV34

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tyrosine-kinase-with-immunoglobulin-like-and-egf-like-domains-2-tie-2/203377>