

Human Plasma Serine Protease Inhibitor (SERPINA5) ELISA Kit | SERPINA5 elisa kit

Human Plasma Serine Protease Inhibitor (SERPINA5) ELISA Kit

CATALOG #
AAA203392

REACTIVITY
Human


PRODUCT SYNONYMS

Plasma Serine Protease Inhibitor (SERPINA5); Human Plasma Serine Protease Inhibitor (SERPINA5) ELISA Kit; SERPINA5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	406
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SERPINA5 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SERPINA5 in the sample, This SERPINA5 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SERPINA5 concentration. The concentration of SERPINA5 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	194018472
NCBI GeneID	5104
NCBI GenBank Nucleotide #	NP_000615.3
UniProt Accession #	P05154
NCBI Official Full Name	plasma serine protease inhibitor preproprotein
NCBI Official Synonym Full Names	serpin family A member 5
NCBI Official Symbol	SERPINA5
NCBI Official Synonym Symbols	PCI; PAI3; PAI-3; PCI-B; PROC1; PLANH3

NCBI Protein Information	plasma serine protease inhibitor
UniProt Protein Name	Plasma serine protease inhibitor
UniProt Gene Name	SERPINA5
UniProt Synonym Gene Names	PCI; PLANH3; PROCI; PAI-3; PAI3; PCI
UniProt Entry Name	IPSP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/plasma-serine-protease-inhibitor-serpina5/203392>