

Human Tumor Necrosis Factor Receptor-Associated Factor 3 (TRAF3) ELISA Kit | TRAF3 elisa kit

Human Tumor Necrosis Factor Receptor-Associated Factor 3 (TRAF3) ELISA Kit

CATALOG #
AAA203579

REACTIVITY
Human


PRODUCT SYNONYMS

Tumor Necrosis Factor Receptor-Associated Factor 3 (TRAF3); Human Tumor Necrosis Factor Receptor-Associated Factor 3 (TRAF3) ELISA Kit; TRAF3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	0.5ng/ml-8ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TRAF3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TRAF3 in the sample, This TRAF3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TRAF3 concentration. The concentration of TRAF3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/tumor-necrosis-factor-receptor-associated-factor-3-traf3/203579>