

## Human Cathepsin L (CATH-L) ELISA Kit | CATH-L elisa kit

Human Cathepsin L (CATH-L) ELISA Kit

CATALOG #  
**AAA203621**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Cathepsin L (CATH-L); Human Cathepsin L (CATH-L) ELISA Kit; CATH-L elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Human                                                             |
| Sequence Length | 306                                                               |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 200 pg/ml-3200 pg/ml                                              |
| Sensitivity     | 10 pg/ml                                                          |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

Related Product Information for CATH-L elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Cath-L in the sample, This Cath-L Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Cath-L concentration. The concentration of Cath-L in the samples is then determined by comparing the O.D. of the samples to the standard curve.

### NCBI AND UNIPROT INFORMATION

|                         |                      |
|-------------------------|----------------------|
| NCBI GI #               | 32394730             |
| NCBI Official Full Name | cathepsin L, partial |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cathepsin-l-cath-l/203621>