

Human N Terminal Probrain B Type Natriuretic Peptides ELISA Kit

Human N Terminal Probrain B Type Natriuretic Peptides ELISA Kit

CATALOG #
AAA203734

REACTIVITY
Human


PRODUCT SYNONYMS

N Terminal Probrain B Type Natriuretic Peptides; Human N Terminal Probrain B Type Natriuretic Peptides ELISA Kit; N Terminal Probrain B Type Natriuretic Peptides elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sensitivity	1 pg/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for N Terminal Probrain B Type Natriuretic Peptides elisa kit

Intended Uses: This N TERMINAL PROBRAIN B TYPE NATRIURETIC PEPTIDES ELISA kit is intended Laboratory for Research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of N TERMINAL PROBRAIN B TYPE NATRIURETIC PEPTIDES in the sample, this N TERMINAL PROBRAIN B TYPE NATRIURETIC PEPTIDES ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus N TERMINAL PROBRAIN B TYPE NATRIURETIC PEPTIDES concentration. The concentration of N TERMINAL PROBRAIN B TYPE NATRIURETIC PEPTIDES in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.