

Human Peroxynitrite Anion ELISA Kit

Human Peroxynitrite Anion ELISA Kit

CATALOG #
AAA203751

REACTIVITY
Human


PRODUCT SYNONYMS

Peroxynitrite Anion; Human Peroxynitrite Anion ELISA Kit; Peroxynitrite Anion elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sensitivity	1 ng/L
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Peroxynitrite Anion elisa kit

Intended Uses: This PEROXYNITRITE ANION ELISA kit is intended Laboratory for Research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PEROXYNITRITE ANION in the sample, this PEROXYNITRITE ANION ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PEROXYNITRITE ANION concentration. The concentration of PEROXYNITRITE ANION in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.