

Human Clostridium Difficile Toxin B ELISA Kit | toxB elisa kit

Human Clostridium Difficile Toxin B ELISA Kit

CATALOG #
AAA203781

REACTIVITY
Human


PRODUCT SYNONYMS

Clostridium Difficile Toxin B; Human Clostridium Difficile Toxin B ELISA Kit; toxB elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	2366
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	25 pg/ml-400 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for toxB elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CLOSTRIDIUM DIFFICILE TOXIN B in the sample, This CLOSTRIDIUM DIFFICILE TOXIN B ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CLOSTRIDIUM DIFFICILE TOXIN B concentration. The concentration of CLOSTRIDIUM DIFFICILE TOXIN B in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	136015
UniProt Accession #	P18177

NCBI Official Full Name	Toxin B
UniProt Protein Name	Toxin B
UniProt Gene Name	toxB
UniProt Synonym Gene Names	tcdB
UniProt Entry Name	TOXB_PEPDI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/clostridium-difficile-toxin-b/203781>