

Human N-Glycolylneuraminic Acid ELISA Kit

Human N-Glycolylneuraminic Acid ELISA Kit

CATALOG #
AAA203961

REACTIVITY
Human

PRODUCT SYNONYMS

N-Glycolylneuraminic Acid; Human N-Glycolylneuraminic Acid ELISA Kit; N-Glycolylneuraminic Acid elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	62.5 nmol/L-1000 nmol/L
Sensitivity	1 nmol/L
CAS	1113-83-3

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for N-Glycolylneuraminic Acid elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of N-GLYCOLYLNEURAMINIC ACID in the sample, This N-GLYCOLYLNEURAMINIC ACID ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus N-GLYCOLYLNEURAMINIC ACID concentration. The concentration of N-GLYCOLYLNEURAMINIC ACID in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.