

Human A Disintegrin And Metalloprotease 7 ELISA Kit | ADAMTS7 elisa kit

Human A Disintegrin And Metalloprotease 7 ELISA Kit

CATALOG #
AAA204111

REACTIVITY
Human


PRODUCT SYNONYMS

A Disintegrin And Metalloprotease 7; Human A Disintegrin And Metalloprotease 7 ELISA Kit; ADAMTS7 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1686
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	12.5ng/ml-200ng/ml
Sensitivity	1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADAMTS7 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of A DISINTEGRIN AND METALLOPROTEASE 7 in the sample, This A DISINTEGRIN AND METALLOPROTEASE 7 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus A DISINTEGRIN AND METALLOPROTEASE 7 concentration. The concentration of A DISINTEGRIN AND METALLOPROTEASE 7 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	38683827
NCBI GeneID	11173

NCBI GenBank Nucleotide #	NP_055087.2
UniProt Accession #	Q9UKP4
NCBI Official Full Name	A disintegrin and metalloproteinase with thrombospondin motifs 7 preproprotein
NCBI Official Synonym Full Names	ADAM metalloproteinase with thrombospondin type 1 motif 7
NCBI Official Symbol	ADAMTS7
NCBI Official Synonym Symbols	ADAM-TS7; ADAMTS-7; ADAM-TS 7
NCBI Protein Information	A disintegrin and metalloproteinase with thrombospondin motifs 7
UniProt Protein Name	A disintegrin and metalloproteinase with thrombospondin motifs 7
UniProt Gene Name	ADAMTS7
UniProt Synonym Gene Names	ADAM-TS 7; ADAM-TS7; ADAMTS-7
UniProt Entry Name	ATS7_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/a-disintegrin-and-metalloprotease-7/204111>