

Human Cluster of Differentiation 36 ELISA Kit

Human Cluster of Differentiation 36 ELISA Kit

CATALOG #
AAA204136

REACTIVITY
Human


PRODUCT SYNONYMS

Cluster of Differentiation 36; Human Cluster of Differentiation 36 ELISA Kit; Cluster of Differentiation 36 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	18
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	1.5ng/ml-24ng/ml
Sensitivity	0.1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Cluster of Differentiation 36 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CLUSTER OF DIFFERENTIATION 36 in the sample, This CLUSTER OF DIFFERENTIATION 36 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CLUSTER OF DIFFERENTIATION 36 concentration. The concentration of CLUSTER OF DIFFERENTIATION 36 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	1646662883
NCBI Official Full Name	cluster of differentiation 36, partial

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/cluster-of-differentiation-36/204136>