

Human Calcium/Calmodulin-Dependent Protein Kinase Type 1D ELISA Kit | CAMK1D elisa kit

Human Calcium/Calmodulin-Dependent Protein Kinase Type 1D ELISA Kit

CATALOG # AAA204155	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Calcium/Calmodulin-Dependent Protein Kinase Type 1D; Human Calcium/Calmodulin-Dependent Protein Kinase Type 1D ELISA Kit; CAMK1D elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	385
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	156.25 pg/ml-2500 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CAMK1D elisa kit	Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Calcium/Calmodulin-dependent protein kinase type 1D in the sample, This Calcium/Calmodulin-dependent protein kinase type 1D Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Calcium/Calmodulin-dependent protein kinase type 1D concentration. The concentration of Calcium/Calmodulin-dependent protein kinase type 1D in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	23273491
-----------	----------

NCBI GeneID	57118
NCBI Official Full Name	Calcium/calmodulin-dependent protein kinase ID
NCBI Official Synonym Full Names	calcium/calmodulin dependent protein kinase ID
NCBI Official Symbol	CAMK1D
NCBI Official Synonym Symbols	CKLiK; CaM-K1; CaMKID
NCBI Protein Information	calcium/calmodulin-dependent protein kinase type 1D

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/calcium-calmodulin-dependent-protein-kinase-type-1d/204155>