

Human Salusin-Beta ELISA Kit

Human Salusin-Beta ELISA Kit

CATALOG #
AAA204337

REACTIVITY
Human


PRODUCT SYNONYMS

Salusin-Beta; Human Salusin-Beta ELISA Kit; Salusin-Beta elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	2.5 pg/ml-40 pg/ml
Sensitivity	0.1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Salusin-Beta elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Human Salusin-Beta in the sample, This Human Salusin-Beta Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Human Salusin-Beta concentration. The concentration of Human Salusin-Beta in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.